

Template DNA: previously genotyped DNA prepared by proteinase K and phenol/chloroform extraction was used at 50ng per 50µl reaction. Genotypes were typically one homozygous A/A, one homozygous C/C and one heterozygote (A/C).

Buffer (1x): 10 mM Tris-HCl (pH 8.3), 1.2 mM or 3.5 mM MgCl₂, 50 mM KCl, dNTPs (each at 100 µM), gelatin at 0.01% (w/v).

Enzyme: AmpliTaq Gold (Perkin-Elmer/ABI) was included in the reaction mix at 2units/50µl reaction.

Page 22, first full paragraph, replace with the following:

Examples 7 and 8

Random coil embodiment and bimolecular embodiment

Scorpion B2731:

fam-AGGTAGTGCAGAGAGTG-mr-h-GAGCCTCAACATCCTGCTCCCCTCCTACTAC (SEQ ID NO: 4); (SEQ ID NO:5)

Scorpion B4249 (no quencher on same molecule)

fam-AGGTAGTGCAGAGAGTG-h-GAGCCTCAACATCCTGCTCCCCTCCTACTAC (SEQ ID NO:6); (SEQ ID NO: 5)

Quencher oligonucleotide (complement of the tail of B4249):

CACTCTCTGCACTACCT-mr (SEQ ID NO:7)

ARMS primer R284-97: TTCGGGGCTCCACACGGCGACTCTCAAC (SEQ ID NO:8)

ARMS primer R283-97: TTCGGGGCTCCACACGGCGACTCTCAAG (SEQ ID NO:9)

Target is the H63D polymorphism of the human hereditary haemochromatosis gene (H/H), B2731 and B4249 are "common" primers to oppose the ARMS primers R283-97, R283-97. Cycling conditions and reaction composition as above. Primers (including *Scorpion* primers) were used at 500nM concentration.

Page 23, first full paragraph, replace with the following:

Example 9

No quencher embodiment

Scorpion B4249 (no quencher)